

# UREA U.V.

## Enzymatic and colourimetric determination of Urea in biological liquids

### TEST SUMMARY

The Urea is being hydrolyzed by means of Urease to Ammonia, which reacts in presence of Glutamate dehydrogenase, with 2-oxoglutarate and NADH to form glutamate and NAD<sup>+</sup>.  
The decrease of absorbance is measured at 340 nm.

### SAMPLES

Serum, plasma (avoid ammonium heparinate).  
Urine diluted 1:100 with deionized water.  
Urea is stable 3 days at 2-8°C.

### REAGENTS

Reagent A: Tris pH 7.60 115 mM, ADP 1.2 mM, NADH 0.25 mM, Urease ≥ 8 KU/lt, GLDH ≥ 800U/lt, stabilizers.

Reagent B: 2-oxoglutarate Sodium 7.5 mM, stabilizers.

Standard: Urea 50 mg/dl; stabilizers and preservatives.

### MATERIALS REQUIRED BUT NOT SUPPLIED

Current laboratory instrumentation. Spectrophotometer UV/VIS with thermostatic cuvette holder. Automatic micropipettes. Glass or high quality polystyrene cuvettes. Saline solutions.

### PRECAUTIONS

Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow. Perform the test according to the general "Good Laboratory Practice" (GMP) guidelines.

### REAGENTS PREPARATION

Add 20 ml of reagent B to a bottle of reagent A.  
Stability: 60 days at 2-8°C.  
Keep at room temperature before use.

### TEST PROCEDURE (STARTER SAMPLE)

Kind of analysis: Fixed time (decreasing)  
Reading time: 30, 90 seconds  
Delay: 30 seconds  
Wavelength: 340 nm  
Temperature: 37°C  
Lightpath: 1 cm  
Zero: Blank reagent

REAGENTS	BLANK	STANDARD	SAMPLE
Working reagent	1 ml	1 ml	1 ml
Incubate at 37 °C for 5 minutes			
Distilled water	10 µl	--	--
Standard	--	10 µl	--
Sample	--	--	10 µl

Record the absorbance after 30 seconds (A1) at 340 nm, record again after exactly 60 seconds (A2).

### CALCULATION

**Serum/Plasma** Urea (mg/dl)  

$$\frac{[A_2 - A_1 \text{ (sample)}]}{[A_2 - A_1 \text{ (standard)}]} \times 50$$
 (standard value)

**Random urine** Urea (mg/dl)  

$$\frac{[A_2 - A_1 \text{ (sample)}]}{[A_2 - A_1 \text{ (standard)}]} \times 50 \times 100$$
 (standard value and dilution)

**24 hours urine** Urea (g/24h)  

$$\frac{[A_2 - A_1 \text{ (sample)}]}{[A_2 - A_1 \text{ (standard)}]} \times 50 \times 100 \times \text{dl urine}$$
 1000

(standard value, dilution factor and diuresis in decilitres)

### EXPECTED VALUES

Adults 10 – 50 mg/dl (1.7 – 8.3 mmol/l)  
 Urine 20 – 35 g/24 h (332 – 580 mmol/24 h)

Each laboratory should establish appropriate reference intervals related to its population.

### NOTE

- If the results are incompatible with clinical presentation, they have to be evaluated within a total clinical study.
- Only for IVD use.

### CALIBRATION/QUALITY CONTROL

It is suggested to perform an internal quality control. For this purpose the following control sera on human base are available on request:

**QN 0050 CH** 10 x 5 ml  
 Control Sera normal values

**QP 0050 CH** 10 x 5 ml  
 Control Sera pathological values

### TEST PERFORMANCE

#### Precision

Intra-assay (n = 30)	Mean (mg/dl)	SD (mg/dl)	CV%
Sample 1	41.23	0.7738	1.88
Sample 2	129.7	0.7497	0.58

Inter-assay (n = 30)	Mean (mg/dl)	SD (mg/dl)	CV%
Sample 1	40.56	1.3308	3.28
Sample 2	130.30	1.2077	0.93

#### Sensitivity/Limit of detection

The method is able to discriminate until 1 mg/dl.

#### Linearity

The method is linear up to 300 mg/dl. If the values are exceeded, it is suggested to dilute the sample 1+9 with saline and to repeat the test, multiplying the results by 10.

#### Methods comparison

A comparison with a commercial available product gave the following results in a comparison on 100 samples:

Urea LTA = x  
 Urea competitor = y  
 n = 31

$y = 0,98577x + 0.38516 \text{ mg/dl}$      $r = 0,9795$

#### Interferences

No interference was observed by the presence of:  
 haemoglobin ≤ 500 mg/dl  
 bilirubin ≤ 28 mg/dl  
 lipids ≤ 600 mg/dl

#### WASTE DISPOSAL

Product is intended for professional laboratories. Waste products must be handled as per relevant security cards and local regulations.

#### PACKAGING

**CODE CC02301 (400 TESTS)**  
 Reagent A 4 x 80 ml (liquid)  
 Reagent B 2 x 40 ml (liquid)  
 Standard 1 x 5 ml (liquid)









### REFERENCES

Falke, H.N. Schubert, G.E. Klin. Wschr.42 (1965).  
 Tietz Textbook of Clinical Chemistry, Second Edition, Burtis-Ashood (1994).  
 HU Bergmeyer – Methods of enzymatic analysis, (1987).

### MANUFACTURER

LTA s.r.l.  
 Via Milano 15/F  
 20060 Bussero (Milan) ITALY  
 Tel: ++39 02 95409034  
 Fax: ++39 02 95334185  
 e-mail: info@ltaonline.it  
 Website: http://www.ltaonline.it

### SYMBOLS

-  IVD Only for IVD use
-  LOT Lot of manufacturing
-  REF Code number
-  Storage temperature interval
-  Expiration date (year, month)
-  Warning, read enclosed documents
-  Read the directions
-  Biological risk

Mod. 01.06 (ver. 4.4 - 12/02/2009)

