

CONTROL SET COPPER / ZINC

2 Levels control solutions for determination of Copper and Zinc in biological fluids

INTENDED USE

Concentrations and activities of components are studied for an excellent use even in automatic system of clinical chemistry, even in manual procedures.

PRESERVATION AND STABILITY

Preserve the product at 2-8°C, stable until expiration date on the package.
Once opened controls are stable for 7 days at 2-8°C.

MATERIALS REQUIRED BUT NOT SUPPLIED

Normal laboratory equipment.

REAGENTS PREPARATIONS

Controls are ready to use.
Bring controls at room temperature before use and verify that there isn't crystalline precipitate on the bottom of the bottle. In precipitate presence shake and/or warm a little bit to re-dissolve.

PRECAUTIONS

- The product contains very dilute acids.
- Avoid contact with skin and mucous membrane.
- Use the normal precautions for lab behaviour.
- For in vitro use only.

WASTE DISPOSAL

Product is intended for professional laboratories. Waste products must be handled as per relevant security cards and local regulations.

PACKAGING

CODE CC02780


Control Level 1 2 x 5 ml (liquid)
Control Level 2 2 x 5 ml (liquid)

REFERENCES



Archives data.

CALIBRATION VALUES

LOT A352.24*

 2028/10






This Table of valueset for lot is applicable to sublots identified by the sequential alphabetical letter (for example A, B, C, etc.) following the lot number.

COMPONENT / METHOD		VALUE	RANGE	UNIT
CONTROL LEVEL 1		Colorimetric tests for the determination of Copper and Zinc in biological fluids		
LOT H290.24A  2028/10	COPPER	113	93 – 134	µg/dl
	ZINC	91.8	75.3 – 108.4	µg/dl
CONTROL LEVEL 2		Colorimetric tests for the determination of Copper and Zinc in biological fluids		
LOT I290.24A  2028/10	COPPER	330	271 – 389	µg/dl
	ZINC	55.4	45.4 – 65.4	µg/dl

MANUFACTURER

LTA s.r.l.
Via Milano 15/F
20060 Bussero (Milan) - ITALY
Tel: +39 02 95409034
Fax: +39 02 95334185
e-mail: info@LTAonline.it
Website: http://www.LTAonline.it

SYMBOLS

- IVD** Only for IVD use
- LOT** Lot of manufacturing
- REF** Code number
-  Storage temperature interval
-  Expiration date (year – month)
-  Warning, read enclosed documents
-  Read the directions
-  Biological risk

Mod. 01.06 (ver. 1.0 – 07/03/2017)

