# PATHOLOGICAL CONTROL SERUM

## SIERO DI CONTROLLO UNIVERSALE

### MATRICE UMANA

**REF** CC03200 10x5 ml

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**REF** CC03200 10x5 ml

#### Intended use

Pathological Control Serum is for use in the quality control of general clinical chemistry methods for the quantitative determination of substrates, electrolytes, lipids, enzymes and proteins. The control is used to monitor accuracy or precision both for manual techniques and assays on automated clinical chemistry analyzers.

#### Summary

Pathological Control Serum is a lyophilized control based on human serum with concentrations and activities in pathological range.

#### Kit components

Lyophilized control serum 10 x 5 ml

#### Reactive lyophilized components

Human serum with chemical additives and tissue extracts of human and animal origin. The origin of biological additives is as follow:

Alkaline phosphatase	calf intestines
ALT/GPT	porcine heart
ASTGOT	porcine heart
Albumin	human plasma
Amylase, total	human saliva
Amylase, pancreatic	porcine pancreas
Cholesterol	human plasma
Cholinesterase	human plasma
Creatine kinase	human, recombinant
γ-GT	porcine kidney
LDH	porcine heart
Lipase	human, recombinant
GLDH	bovine liver
Total protein	human plasma

The concentrations and activities are lot specific. The exact values and ranges are listed in the enclosed values sheet.

#### Precautions and warnings

For in vitro diagnostic use.

All human material should be considered potentially infectious. All products derived from human blood are prepared exclusively from the blood of donors tested individually and shown by FDA approved method to be free from HBsAg and antibodies HCV and HIV.

However, as no testing method can rule out the potential risk of infection with absolute certainty, the material should be handled with the same level of care as a patient specimen. In the event of exposure, the directives of the responsible health authorities should be followed.<sup>1,2</sup>

#### Disposal

The product is for use within professional laboratories.

Disposal of all waste material should be in accordance with local guidelines. P501: Dispose of contents according to national/international regulations.

#### Handling

Carefully open one bottle of serum, avoiding the loss of lyophilizate, and pipette exactly 5.0 ml of distilled / deionized water. Carefully close the bottle and dissolve the contents completely by occasional gentle swirling within 30 minutes. Avoid the formation of foam.

#### Storage and stability

Store at 2-8°C.

Criterion of the stability data specified by the manufacturer: recovery within ± 13% of initial value.

Stability of lyophilized control serum at 2-8°C up to the stated expiration date. Stability in reconstituted control serum at 2-8°C: 5 days. Stability in reconstituted control serum when frozen once: 28 days (-20°C ± 5°C).

The possible appearance of a slight green coloration has no effect on the recovery of the values. Store control tightly capped and protected from light when not in use.

#### Test procedure

Use supplied materials as provided. Additional materials required: reagents and clinical chemistry analyzers, General laboratory equipment, distilled/deionized water.

#### Assay

Dispense the required volume into a sample cup and analyze in the same way as for the samples. The control should be run daily in parallel with the patient samples and after every calibration. The control intervals should be adapted to each laboratory's individual requirements. Each laboratory should establish QC procedures that conform with local, state and/or federal regulations or accreditation required.

#### Assigned values and ranges

The target values were determined using the method stated in the electronically available or enclosed value sheets. The target value specified is the mean of all values obtained. The corresponding control range is calculated as the target value ± 3 SD. Determinations were made using the testing reagents and analysis systems available at the time.

Results should be within the defined ranges.

Each laboratory should establish guidelines for corrective measures to be taken if values fall outside the range.


#### References

- Occupational Safety and Health Standards: bloodborne pathogens. (29 CFR Part 1910.1030). Federal Register.
- Directive 2000/54/EC of the European Parliament and Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work.

#### Manufacturer

LTA srl, via Milano 15/F, 20060 Bussero (MI) – Italia

#### Symbols

**IVD** Only for IVD use, **REF** Code number, **LOT** Lot of manufacturing,  Expiration date,  Storage temperature interval,  Warning, read enclosed documents,  Read the directions,  Biological risk

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