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**REF** CC03250 10x3 ml

#### Intended use / Summary

Calibrator Serum is a human based multiparametric calibration serum. The concentration and activities have been selected to ensure optimum calibration of automatic analyzers.

#### Kit components

Calibrator Lyophilized 10 x 3 ml

#### Reactive lyophilized components

Human serum with chemical additives and tissue extracts of tissue of human and animal origin. The origin of biological additives is as follow:

ALT/GPT	porcine heart
AST/GOT	porcine heart
Albumin	bovine plasma
Amylase, total	porcine pancreas
Amylase, pancreatic	porcine pancreas
Cholesterol	bovine plasma
Cholinesterase	human serum
Creatine kinase	rabbit muscle
γ-GT	porcine kidney
GLDH	bacterial, recombinant
LDH	porcine heart
Lipase	human pancreas (recombinant)
Triglycerides	chicken egg yolk

#### Calibration values

The calibration values were determined using the method mentioned in the enclosed value sheet. Determinations were performed under strictly standardized conditions on Roche analyzers using Roche and LTA reagents and Roche master calibrators. The calibration values were obtained via single assays performed in different laboratories in several independent series. The calibrator value specified is the mean of the values obtained.

#### Precautions and warnings

For in vitro diagnostic use.

All human material should be considered potentially infectious. All products derived from human blood are prepared exclusively from the blood of donors tested individually and shown to be free from HBsAg and antibodies to HCV and HIV by FDA approved methods.

However, as no testing method can rule out the potential risk of infection with absolute certainty, the material should be handled with the same level of care as a patient specimen. In the event of exposure, the directives of the responsible health authorities should be followed.<sup>(1,2)</sup>

#### Disposal

The product is for use within professional laboratories.

Disposal of all waste material should be in accordance with local guidelines. P501: Dispose of contents according to national/international regulations.

#### Handling

Carefully open one bottle of calibrator, avoiding the loss of lyophilizate, and pipette exactly 3.0 ml of distilled / deionized water. Carefully close the bottle and dissolve the contents completely by occasional gentle swirling within 30 minutes. Avoid the formation of foam.

#### Storage and stability

Store at 2-8°C. Stability of lyophilized calibrator at 2-8°C up to the stated expiration date. Criterion of the stability data specified by the manufacturer: recovery within ± 5% of initial value.

Stability of components in reconstituted calibrator:

at 15-25 °C	8 hours
at 2-8 °C	2 days
at (-15)-(-25) °C	28 days (when frozen once)

#### Exceptions:

Stability of total bilirubin in reconstituted calibrator (stored protected from light):

at 15-25 °C	6 hours
at 2-8 °C	24 hours
at (-15)-(-25) °C	14 days (when frozen once)

Stability of direct bilirubin in reconstituted calibrator (stored protected from light):

at 15-25 °C	3 hours
at 2-8 °C	8 hours
at (-15)-(-25) °C	14 days (when frozen once)

Store calibrator tightly capped and protected from light when not in use.

#### Test procedure

Use supplied materials as provided. Additional materials required: automatic micropipette, reagents and clinical chemistry analyzers, General laboratory equipment, distilled/deionized water.

#### Assay

Dispense the required volume into a sample cup and use as advised.

#### References

- Occupational Safety and Health Standards: bloodborne pathogens. (29 CFR Part 1910.1030). Federal Register.
- Council Directive (2000/54/EC). Official Journal of the European Communities. No. L262 from Oct. 17, 2000.

#### Manufacturer

LTA srl, via Milano 15/F, 20041 Bussero (MI) – Italia

#### Symbols

**IVD** Only for IVD use, **REF** Code number, **LOT** Lot of manufacturing,  Expiration date,  Storage temperature interval,  Warning, read enclosed documents,  Read the directions,  Biological risk

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