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BIOCHEMISTRY CONTROL URINE

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Controllo / Control 1 x 5 ml (liofilo)

INTENDED USE

This product is intended for in vitro diagnostic use in quality control on clinical chemistry systems.

Urices are used to check accuracy and / or precision both in automatic analyzers and in manual use.

Control urine provided with suggested ranges of acceptable values.

COMPOSITION

Human urine lyophilized without preservatives which might interfere with the tests.

STORAGE AND STABILITY

Store the non-reconstituted product at 2-8°C, stable until the expiration date indicated on the package.

The reconstituted calibrator is stable 7 days at 2-8°C.

For longer periods, aliquot the control and freeze it. Avoid repeated freezing.

MATERIALS REQUIRED NOT PROVIDED

5 ml micropipette. Distilled water

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Open the vial very carefully, avoiding any loss of the lyophilized material.

Pipette exactly 5 mL of distilled water.

Close the vial with and let stand for 20 minutes at room temperature.

Swirl gently, avoiding the formation of foam, to ensure complete dissolution of contents.

Use the reconstituted control as indicated in the instructions for the use of the analysis to be performed.

ASSIGNED VALUES

The concentrations and activities of the components must refer to each individual lot. The exact theoretical values, the intervals and the traceability of the individual components are contained in the attached sheets.

The suggested intervals of acceptable values have been calculated from previous experience in interlaboratory variability and are given for orientation only. Each laboratory should establish its own precision parameters.

PRECAUTIONS

• *All human material must be considered potentially infectious. All products of human origin are prepared exclusively from donors individually tested for HBsAg, anti-HCV and HIV antibodies and NON-REACTIVE results.*

Since it is not possible to safely exclude the risk of infection, it is necessary to handle the material with the same precautions adopted for the samples taken by the patients. In the case of an exposure, proceed according to the specific health indications.

• *For in vitro diagnostic use only.*

WASTE DISPOSAL

Product is intended for professional laboratories. Waste products must be handled as per relevant security cards and local regulations.

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LTA s.r.l.
Via Milano 15/F
20060 Bussero (Milano) ITALY
tel. +39 02 95409034 - fax. +39 02 95334185
e-mail. info@LTAonline.it - website. http://www.LTAonline.it

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The enzymes values (U/L e μ kat/L) are for a incubation temperature of 37°C.

C-RSE/IFF: Traceable to the reference system as described by the IFCC Committee on Reference Systems for Enzymes.