# Bence Jones Ouchterlony

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# Bence Jones Ouchterlony

## Bence Jones's proteinuria determination by double radial immunodiffusion according to Ouchterlony technique

### INTENDED USE

The device consists of a plate with agarose gel, in the holes of which the antibody and antigen are deposited, respectively; The latter is qualitatively detected through the formation of an immunoprecipitation band in the gel.

The test determines the total (free and bound) kappa and lambda light chains in urine samples. The device is intended for professional use only. For in vitro diagnostic use.

### TEST SUMMARY

The antibody and the antigen dispensed into two different wells diffusing through the gel create a concentration gradient. Between the two wells, where the antibody and antigen concentration are equal there is a formation of immunoprecipitation band.

The test determines the total light chains (free and bound).

### SAMPLES

Urine. Stability 7 days a 2-8°C, 2 months at -18°C.

### REAGENTS

#### Plates

Agarose gel 0,9%, Tris Buffer 50 mM pH 7.4, colouring, PEG 6000 3%.

#### Kappa Antiserum

Goat anti human kappa chains free and bound serum, Tris Buffer 50 mM pH 7.4 preservative and stabilisers.

#### Hazard statements (H):

EUH032 – Contact with acids liberates very toxic gas

#### Lambda Antiserum

Goat anti human Lambda chains free and bound serum, Tris Buffer 50 mM pH 7.4 preservative and stabilisers.

#### Hazard statements (H):

EUH032 – Contact with acids liberates very toxic gas

#### Positive Control

Stabilised human solution with Kappa chains Lambda chains free and bound.

#### Hazard statements (H):

EUH032 – Contact with acids liberates very toxic gas

### REAGENTS PREPARATION AND STORAGE

All reagents are ready to use.

The reagents are stable until expiration date on the label at 2-8°C.

### MATERIALS REQUIRED BUT NOT SUPPLIED

Micropipette to 10-40 µL, incubator 37°C, normal laboratory instrumentation.

### PRECAUTIONS

Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow. Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

### PROCEDURE

Allow all reagent and samples to come to room temperature. And mix well by gentle inversion before use. All unused reagents should be refrigerated as soon as possible after use.


Open the plate and eventually wait the evaporation of moisture present in the wells.

Properly identify the plate and lateral well (kappa well or lambda well).

Pipet 40 µL of urine into central well and 10 µL of Kappa antibody in lateral well and 10 µL Lambda antibody in other well according to the marker on the bottom.

Incubate in wet chamber at 37°C.

If the sample's concentration is high, it's possible to read the results after 6-9 hours of incubation; otherwise, results must be read after 9-15 hours.

### RESULTS INTERPRETATION

The absence of any immunoprecipitation band means that the sample is negative for Bence Jones.

The presence of immunoprecipitation between the sample well and the well of the Kappa antiserum indicates the presence in the sample of Kappa chains.

The presence of immunoprecipitation between the sample well and the well of the Lambda antiserum indicates the presence in the sample of Lambda chains.

Both Kappa chains and Lambda chains detected by test, may be polyclonal, therefore give a negative result at immunofixation test.

A semiquantitative results can be calculated looking at the position of the band; band will be close to the sample well if the concentration is low otherwise it'll be close to the antibody well.

If the test is performed like a preliminary screening to contrastographic assay it is generally accepted a cut off of 50 mg/dL.

However, each laboratory has to determinate its cut off value.

### NOTES

• The Ouchterlony Test should be used for screening purposes as, in patients with moderate to high proteinuria, the test may detect false positives; To confirm the diagnosis, a further test with a different method is necessary.

• As with any diagnostic procedure, if the results are inconsistent with the clinical presentation, the physician should evaluate the data obtained using this test in light of other clinical information.

### WASTE DISPOSAL

This product is made to be used in professional laboratories. Please consult local regulations for a correct waste disposal.

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REF	UK01001	(10 TESTS)
	Piastrine per doppia immunodiffusione / Double immunodiffusion plates 10 pcs	
	Antisiero Kappa / Kappa Antiserum	1 x 300 µL
	Antisiero Lambda / Lambda Antiserum	1 x 300 µL
	Controllo positivo / Positive Control	1 x 100 µL

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**NOTE: MODIFICATIONS HIGHLIGHTED WITH GRAY BACKGROUND**

### CALIBRATION

The Kit is supplied with a Positive Control.

### TEST PERFORMANCE

#### Precision

Intra assay (n= 60)	+	-
Bence Jones LTA	20	40
Test competitor	20	40

#### Sensitivity limits

The sensitivity is 5 mg/dL for Kappa chains and 5 mg/dl for Lambda chains.

### EXPECTED VALUES

The analytes should be absent or present at a concentration lower than cut off of analysis.

**LTA Sri**

Via Milano 15/F – 20041 Bussero (Milano) Italia  
Tel. +39 02 95409034 – Fax +39 02 95334185  
e-mail: info@LTAonline.it – web: www.LTAonline.it