Faecal Parasites Concentrator

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
Fig. 1:

Faecal Parasites Concentrator

Method for concentration of helminth larvae and eggs, protozoan cysts and coccidian oocysts for microscopic examination

PRINCIPLE

The methodology is not an easy alteration of the traditional methods as Formalin Ether by Ritchie, but it is an easy and effective system that allows the immediate separation of parasites from faeces, without the slightest damage for them.

The system concern in only one test-tube which contains substances with high molecular weight, mixed with solutions that operate on faecal material melting it and, in case of presence, making easier the release of eggs, larvae and oocysts. Furthermore, there is a colouring matter in every test-tube that allows a better underlining of eggs, larvae, nucleus, cytoplasm, membrane, and cilia without alter them.

INTENDED USE

LTA faecal parasites concentrator is a method intended for the detection of parasites in the faeces. This reagent is used to concentrate helminth larvae and eggs, protozoan cysts and to recover oocysts like *Giardia lamblia*, *Entamoeba*, *Cryptosporidium parvum* and *Isoospora belli* etc. from stool samples before microscopic analysis. The device is intended for professional use only. For in vitro diagnostic use.

REAGENTS

Concentrator test-tubes: cromotrope 2R, saline solution.

Hazard statements (H):

H302 – Contact with acids liberates very toxic gas.

H410 – Safety data sheet available on request.

MATERIALS REQUIRED NOT SUPPLIED

Normal laboratory equipment. Vortex, centrifuge, microscope.

PRECAUTIONS

Reagent may contain not reactive and conservative components. It is opportune to avoid contacts with skin and ingestion. Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

SAMPLES

Fresh faeces. Stability: 24 hours at 2-8°C.

REAGENTS PREPARATION

Reagent is ready to use. Store at 4-30°C until expiration date on the label.

PROCEDURE

Take one Concentrator test-tube, remove the screw cap gifted with a dosimeter small shovel, melt 2 spoonfuls of faeces in the test-tube's liquid. Mix on vortex and centrifuge for 5 minutes at low turn number (700/900).

N.B. Never exceed those limits.

After centrifugation, draw 2 drops of supernatant from the liquid's surface and, 2 drops from the bottom with a Pasteur pipette. Then proceed with the sample's microscope survey, placing the sample to examine on a slide recovered with another slide.

It is recommended to observe with particular attention also the slide perimeter. See Fig.1 for the illustrated procedure.

NOTES

- If the results are incompatible with a clinical presentation, the physician should evaluate the data obtained within a total clinical study.

TEST PERFORMANCE

Methods comparison.

A comparison with an available commercial method gave the following results on 65 samples compared:

LTA Faecal Parasites = y
Competitor's Faecal Parasites = x
n = 65

y = x r = 1,0

Limitations of the method

Due to intermittent parasites elimination, negative results should be repeated from new samples collected on different days. It is advisable to examine at least three stool samples collected on different days.

WASTE DISPOSAL

Product is intended for professional laboratories. Waste products must be handled as per relevant security cards and local regulations.

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REF VK03500 (50 TESTS)
Faecal parasites concentrator 50 provette / 50 test tubes

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