

FEBRILE NEGATIVE CONTROL

Proteinus control solution for agglutination reactions with bacterial suspensions

REAGENTS

Negative control: Non-reactive Proteinus bovine solution for agglutination reactions with bacterial suspensions in tube, in microplate and in slide.

REAGENTS PREPARATION

Reagents are ready to use.
Stability: until expiration date on label, stored at 2-8°C.
Do not freeze.

PRECAUTION

Reagent may contain some non-reactive and preservative components.
It is suggested to handle carefully it, avoiding contact with skin and swallow.
Perform the test according to the general "Good Laboratory Practice" (GPL) guidelines.

NOTE

Only for IVD use.

WASTE DISPOSAL

Product is intended for professional laboratories.
Waste products must be handled as per relevant security cards and local regulations.

PACKAGING

CODE BS00020

Febrile Negative control 1 x 1 ml









REFERENCES

Edward J Young. Clinical Infectious Diseases 1995; 21: 283-290.
David A et al. Current Opinion in Infectious Diseases 1994; 7: 616-623.
David R et al Current Opinion in Infectious Diseases 1993; 6: 54-62.
Chart H. Serodiag Immunother Infect Diseases 1995; 7: 34-39.
Chart H. Serodiag Immunother Infect Diseases 1995; 7: 30-33.
Bradley D Jones. Annu Rev Immunol 1996; 14: 533-561.

MANUFACTURER

LTA s.r.l.
Via Milano 15/F 20060
Bussero (Milano) ITALY
Tel: ++39 02 95409034
Fax: ++39 02 95334185
e-mail: info@ltaonline.it
Website: <http://www.ltaonline.it>

SYMBOLS

-  Only for IVD use
-  Lot of manufacturing
-  Code number
-  Storage temperature interval
-  Expiration date (year, month)
-  Warning, read enclosed documents
-  Read the directions
-  Biological risk

Mod. 01.06 (ver 2.1 - 05/12/2005)

